

A Few Calculation Hints (extra sig figs shown)

25.0 mL of a 0.175 M KIO_3 solution is titrated to the end-point with 28.0 mL of a $\text{Na}_2\text{S}_2\text{O}_3$. What is the molar concentration of the $\text{Na}_2\text{S}_2\text{O}_3$?

Remember: 1 mol IO_3^- reacts with 6 mol $\text{S}_2\text{O}_3^{2-}$

Answer: 0.9375 M

25.0 mL of this $\text{Na}_2\text{S}_2\text{O}_3$ solution was used to titrate a solution containing 3.2 mL of the original bleach solution.

| | | |
|----|---|---------|
| 1 | volume of original bleach sample (mL) titrated | 3.2 |
| 2 | mass, in g, of original bleach sample titrated* | 3.47 |
| 3 | final buret reading | |
| 4 | initial buret reading | |
| 5 | mL $\text{Na}_2\text{S}_2\text{O}_3$ added | 25.0 |
| 6 | mol $\text{Na}_2\text{S}_2\text{O}_3$ added | 0.02344 |
| 7 | mol ClO^- reacted | 0.01172 |
| 8 | mol "available chlorine" in titrated sample | 0.01172 |
| 9 | grams "available chlorine" in titrated sample | 0.8309 |
| 10 | percent "available chlorine" in original sample | 23.94 |
| 11 | average percent "available chlorine" in original sample | — |
| 12 | average percent, by mass, of NaOCl in original sample | 25.14 |

*assume the density of the bleach sample is 1.084 g/mL

Hints:

**Remember: available chlorine (Cl_2) relationship: 1 mol Cl_2 produces 1 mol ClO^-

Row 7: ? mol $\text{ClO}^- = 0.02344$ mol $\text{Na}_2\text{S}_2\text{O}_3$

2 mol $\text{S}_2\text{O}_3^{2-}$ reacts with 1 mol I_2 ; 1 mol I_2 is produced when 1 mol ClO^- reacts.

Row 8: ? mol $\text{Cl}_2 = 0.01172$ mol ClO^-

Row 9: ? g $\text{Cl}_2 = 0.01172$ mol Cl_2

Row 10: ? g $\text{Cl}_2 = 100$ g liquid bleach

Row 11: Average of Row 10 results

Row 12 (in lab this uses row 11 average): ? g NaOCl = 100 g liquid bleach